

Nitrous Oxide Conscious Sedation

ClearView SINGLE-USE NASAL HOODS CO₂ Capnography Nasal Hood

Designed for the comfort and safety of your patients and for monitoring the effects of their sedation therapy.

Secure Female Luer Lock Sample Line Connection



Connect a capnograph monitor and verify patient end-tidal CO₂

Single-use nasal hood for the delivery and scavenging of nitrous oxide and for end-tidal CO₂ capnograph monitoring. Designed to be use with a CO₂ capnograph monitor system, ClearView CO₂ allows clinicians to comply with state and national healthcare organizations' guidelines regarding capnograph monitoring during moderate to deep sedation and general anesthesia.

- Female Luer lock allows for easy, secure and leak-free connection to capnograph sample line
- Ergonomic, integrated sample line extension follows the flow of the scavenging circuit's tubing, remaining kink-free and away from the field of view
- Transparent outer hood provides visual of nose breathing – fogging appears in hood upon exhalation, dissipates upon inhalation
- Hood-in-hood design improves nitrous scavenging efficiency
- Low profile improves clinician's view and access to mouth
- Latex-free material protects patient and staff who have latex allergies
- Individually wrapped nasal hood offers patient assurance that they are the only user
- Single-use:
 - Minimizes the risk of cross contamination associated with improper cleaning and sterilization
 - Cost effective – reduces staff time spent on sterilization procedures



Compliance with the 2014 AAOMS guidelines requires the use of CO₂ capnograph monitor when administering moderate to deep sedation and general anesthesia.



Latex-free

PEDO SIZES:

- Birthday Bubblegum 43037-12
- Unscented 43037-9

ADULT SIZES:

- French Vanilla 43035-14
- Unscented 43035-9

12 nasal hoods per box

ACCUTRON Inc.
www.accutron-inc.com

Relentless Innovation. Creative Solutions.

ClearView

SCAVENGING CIRCUITS

Accutron scavenging circuits help to secure a safe and productive work environment for patients, doctor and staff.

- X-ray compatible circuitry saves time if images are required during treatment
- Easy-to-clean, autoclavable circuitry protects against cross-contamination
- Components and rubber goods are latex-free
- Made in the U.S.A. for maximum quality standards
- Bright white finish complements clinical operatory environments
- Flexible system accommodates most brands of flowmeters



Order No. 43003

Check your scavenging circuit's breathing bag and corrugated tubing for tears and holes on a regular basis.

ClearView Standard Scavenging Circuit

43003 ClearView Scavenging Circuit II

Includes scavenging circuit, rubber goods (breathing bag and corrugated tube), and a complementary sampling of ClearView Single-Use Nasal Hoods

43005 ClearView Scavenging Circuit I

Includes scavenging circuit and a complementary sampling of ClearView Single-Use Nasal Hoods

ClearView Scavenging Circuit for Remote Flow System (RFS) and Accu-Vac

43004 ClearView Scavenging Circuit II

Includes scavenging circuit, rubber goods (breathing bag and corrugated tube), and a complementary sampling of ClearView Single-Use Nasal Hoods

43006 ClearView Scavenging Circuit I

Includes scavenging circuit and a complementary sampling of ClearView Single-Use Nasal Hoods

Vacuum Controllers

36855 Remote Flow System S (for Ultra PC % Cabinet Mount Flowmeter)

56855 Remote Flow System D (for Digital Ultra Flushmount Flowmeter)

56955 Accu-Vac for Digital Ultra™ Flowmeter

All Accutron products are available through your authorized dental distributor.



1733 Parkside Lane • Phoenix, AZ 85027
Toll-free: 800.531.2221 • Local Phone: 623.780.2020 • Fax: 623.780.0444

www.accutron-inc.com



Scan QR code

**RFS Chairmount
Kits are available
for the Remote
Flow System**



Forest chair shown